

FDA Perspective: “The Role of Basic Science in the Design of Clinical Trials for CNS Devices”

Victor Krauthamer, Ph.D.

Office of Science and Engineering Labs
Center for Devices and Radiological Health


February 26, 2008

“The opinions and/or conclusion expressed are solely those of the author and in no way imply a policy or position of the Food and Drug Administration.”




Some CNS Devices


Deep Brain Stimulator for Tremor



Spinal Cord Stimulator for Pain



Cochlear Implant for Hearing



Functions of Basic Science and Engineering in CNS Device Trials

- Is a clinical trial necessary? *85 - 90% of clearances for all devices do not need clinical trials*
 - What type of clinical trial is necessary? *e.g., type of control group*
 - Does prior knowledge of mechanism provide confidence?
 - What can be said about safety and efficacy before clinical trial?
 - What are the clinical endpoints for safety and effectiveness?
Causal-path biomarkers → surrogate endpoints
 - Can trials be focused, and less burdensome?
 - What can be predicted beyond the empirical trial? *side effects, interactions with drugs, help interpret results*
-

When Are Clinical Data Needed?

To support:

- PMA, PDP or HDE (almost always)
- 510(k) (10 – 15%*, usually involving new technology that could affect safety of effectiveness or a new population)
- New indication for an *approved* device (e.g., new indications for DBS)
- Significant change to device, especially Class III devices

*for all devices

When is a randomized, double-blinded, placebo-control trial most necessary?

- When no prior information exists – e.g., no previous experiments, when causal path is not known
 - When prior information cannot be trusted – e.g., underpowered or poorly controlled trials, especially when there is a high likelihood of a placebo effect or regression to the mean
-
-

When clinical trials may not be justified?

- When Causal Effects are Known
 - to evaluate a scalpel: cuts because it is sharp
 - to evaluate a sterilizer: temperature and pressure
 - mathematical models used to support safety of a cochlear implant in MR imaging in a PMA
- Codified use of prior information: section 510(k)

The better the prior information, the less justified the randomized placebo-control study.

- Ethical Reasons – Declaration of Helsinki limits placebo controls
-

Fundamental science is provides the necessary context outside of a clinical trial.

Is the empirical trial enough?

- Willingness to infer causality
 - Willingness to accept uncertainty, $\alpha \leq 5\%$
 - Trails are never perfect
 - Real practice is different from clinical trial
 - Demographic differences
 - Duration of study
 - Patient drop-out
 - Blinding issues
 - Unknown covariates
-

concept from Terrence Blaschke, Stanford U.

em·pir·i·cism

Function: *noun*

Date: 1657

1 a : a former school of medical practice founded on experience without the aid of science or theory **b** :

QUACKERY, CHARLATANRY

2 a : the practice of relying on observation and experiment especially in the natural sciences **b** : a tenet arrived at empirically

3 : a theory that all knowledge originates in experience

Simpson's Paradox

| Combined | Recovered | Did not recover | N | Recovery Rate |
|-----------------|-----------|-----------------|----|---------------|
| drug | 20 | 20 | 40 | 50% |
| no-drug | 16 | 24 | 40 | 40% |

| Males | Recovered | Did not recover | N | Recovery Rate |
|--------------|-----------|-----------------|----|---------------|
| drug | 18 | 12 | 30 | 60% |
| no-drug | 7 | 3 | 10 | 70% |

| Females | Recovered | Did not recover | N | Recovery Rate |
|----------------|-----------|-----------------|----|---------------|
| drug | 2 | 8 | 10 | 20% |
| no-drug | 9 | 21 | 30 | 30% |

*Uses of Fundamental Science
within the Clinical Trial*




Preclinical Science - Electrical Stimulation Safety in CNS

- Heating (Energy - Joules)
 - Electrochemical Toxicity (Charge Density, Charge Balance)
 - Dielectric Breakdown (tissue E-Field)
 - Changes to ion channels (E- field)
 - Excitotoxicity (Pulse Repetition Rate, Activation Function ~ Tissue E-field)
 - Physiological Interference (E-field, Repetition Rate)
-



Safety Endpoints

- Standard Ones
 - Neurotoxicity
 - Histopathology (typically staining of cell bodies, glia and myelin)
 - Nonstandard
 - changes to dendrites and axon terminals
 - alteration or loss of normal function
-


Histology and Brain Stimulation, Nissl Stain




Limits of Standard Histology



In vitro study of
stimulation:
Human
Neuroblastoma
Cells Stimulated
for 1 Hour




Ca²⁺ Elevation - Metabolic Overload, 3hrs




25 μ


0.15-0.40



0.20-0.65



0.45-1.40



2.00-4.50



Clinical Safety in Cochlear Implants

- Absence of pain
 - Long-term changes
 - Electrical threshold
 - Comfort levels
 - Dynamic range (maximum comfort level minus threshold)
 - Electrode impedance changes
 - Neural Function
-

Deep Brain Stimulation – 1997 Approval


| | |
|------------------------|-------------------------------|
| Electrode Surface Area | 0.06 cm ² |
| Pulse Duration | 4500 microseconds |
| Amplitude | 10.5 V, or 21.0 mA into 500 Ω |
| Charge per Phase | 9.45 μCoulombs |
| Charge Density | 150 μCoulombs/cm ² |
| Frequency | 185 Hz |

Problems

1. High charge density possible.
2. Small lesions seen in autopsy.
3. Progressive impairment in 10% of patients.

Solutions

1. The sponsor has agreed to provide programmer warning and labeling for parameters which result in charge densities greater than 30 μCoulombs/cm²
2. Post-approval autopsy studies
3. Post-approval effectiveness study



Looking Forward

- CNS causal-path biomarkers needed – correlates are not good enough.
 - Are biomarkers really needed for affective disorders when patients feelings and function are paramount?
 - How can we compare between trials?
 - Can “placebo” effects be used effectively by device treatments?
-